



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT

Operation Table Hydraulic

Item Code: HF26

High quality operating table with built in multiple functions making it suitable for general surgery Functional design & easy to use hydraulic adjustment Wide accessory range for specialized procedures 3 section table top frame of Stainless steel Antistatic PU Integral mattress Height around 730 - 1040 mm. Trendelenburg 25 -- +25 degree Side tilt +/- 15 deg with gas spring Back section -15 - 70deg with gas spring Leg section - 90 -0 deg with gas spring Headrest -40 --+30 deg gas spring assisted Operational area around 500X1950 Height Capable to rise up to 1.2 meter. supplied with accessory like I.V rod, Lithotomy Legs etc.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com