



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Ophthalmic Operating Microscope

Item Code:

Rajas Operating Microscope - Surgical Operating Microscope Cold Light provides high intensity with flexibility in operating. It is focused with built in manual 3 step changer for continuous adjustment. It is based on Optical Light Guide System and checks the interior part of the eye.

Eye Pieces : 12.5x / 10x / 16mm Dioptic Adjustment

Magnification Changes : Built In Manual 3-Step Changer

Total Magnification : 7.5x, 12.5x, 20x

Field Of View : 37.8mm, 23.5mm, 14.5mm

Working Distance : 200mm

Interpupillary Distance: 55 mm to 75 mm

Adjustment

Type : By Optical Light Guide System

Light Source : 24V/250Watts Halogen Lamp

Illumination Control : Continuous Adjustment

Field Of Illumination : 70mm

Filter : Cobalt Blue, Red Free, Heat Absorbing Filter

Type : Floor Stand

Vertical Range Of: 400 mm

Counterbalanced Arm

General Information

Base Size : 600mm x 550mm

Weight : 50Kg.

Height : 1550 mm

Electrical Power : AC 220V, 240V, 50 / 60 Hz

Power Consumption : 220 VA / 110V is also Available

Yes! I am interested

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar

Ambala Cantt- 133004 (Haryana) India

E-mail : rajasmedical@gmail.com

www.rajasindia.com