



## Ophthalmic Refraction Unit

Item Code: OE-05

Need-based, highly economical combination unit Maximum operative functionality; minimum space consumption Comprehensive, easy-to-use console panel control with memo key. 165-degree reclining, gas spring-driven chair. Smooth, Good quality DC motor-driven elevation. Sliding table with space for two instruments. older for hand instrument Dust-free, enclosed Trial Lens Rack with straight, unobtrusive, in & out sliding movement. Sturdy Lamp Pole to support multi-instrument mounting. Cross-integrated Chair Raising System & Sliding table Top. Comfortable Polyurethane seat. Moulded plastic-covered base area. Broad, square format platform for stable mounting & jerk-free performance. Chair Min Height 550 mm Chair Max Height 750 mm Motor Movement 18mm/Sec Instrument Table Height 900mm Instrument Table Length 32" Instrument Table Width 16" Sliding Moment 340mm Extn Table Size 16"x16" Lamp pole height 2210mm Chair Rotation 180° Trail Box Size (526x355x65)mm Width of Chair Unit 1170+ 340 mm

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)