



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Over Bed Table Fix Height

Item Code:

over bed table has been widely appreciated by our clients worldwide. Available in overall approximate size of 1120 mm x 450 mm, W x 965 mm, it consists of mild steel tubular frame work mounted on four twin wheels. salient features are as follow : The diameter of the non rusting castor wheel is 50mm Durable The top board is pre-laminated and are available in the size of 760 mm x 450mm Available in different KDC Completely pretreated and powder coated.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com