



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Over Bed Table

Item Code:

Hospital over bed table is used for patients for taking breakfast and lunch. The laminated board top in two section measures 820x360 mm. The tubular frame work and adjustable height by means of a gas spring is easy to use and provides comfort to the patients and staff for close observation. construction of the table made of M.S. square pipe finished with duly powder coated paint finish unit moves on wheels two with break and two without break.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com