



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Pediatric Bed

Item Code:

We are instrumental in providing our clients a comprehensive range of Pediatric Bed Cot.

- (a) Overall approx Size: 1450mm x 650mm x 600 mm
- (b) Frame work made of rectangular M.S. tubes
- (c) Perforated CRC Sheet duly machine press bent
- (d) M.S. Powder Coated Tubular Head & Foot Bows of equal height
- (e) Drop Down Safety side railings on both sides.
- (f) Bed Mounted on PVC Stumps
- (g) Pre Treated and Epoxy powder Coated.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)