



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Plasma Expressor Electrical

Item Code: BB-14

Microprocessor controlled, designed for automatic separation of blood components from the blood bag. Automatic Control, Audiovisual Alarm, and provided with infra-Red Sensor. Clamping Motor activated, Mild steel cabinet, white powder coated able to withstand corrosion, wear & tear. The Acrylic Compression plate designed to exert uniform pressure on the blood bag. 230/110V AC.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com