



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Platelet Agitator

Item Code: B-B-07

Rajas Platelet agitator specially designed to maintain, Platelets in continues motion at a fixed speed of 70 RPM. Maximum 36 standard Platelet bags can be accommodate at a time outer size 46x39x38 cms construction made of steel with heavy base duly powder coated paint finish, speed regulator on and off swith, rubber stamps to give a stablelity on floor, mesh made of stainless steel with stainless steel side support, 5 to 6 shelves provided to fill the blood platlets for agitation top easy to clean. DC motor for non stop trouble free working on the speed of 70 to 80 RPMs which give smooth movement to agitator.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com