



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Platelet Incubator

Item Code: BB-09

Rajas Plated Incubator provides an ideal strong environment for the Platelets. This unit designed to maintain a stable and inform temperature of 22.0°, the user friendly includes an acrylic inner door, and see through outer door, simple calibration is possible via the control panel. This, coupled with a visual and audible alarm system, enureses that your samples are both safe and secure, but also easy to access when you need them. Other features of this unit include Micro Controlled Temperature controller with Audio Visual Alarm for high & low temperature Display of set and process value. Extra - strength non CFC refrigeration system The unit is mounted on castor wheel for easy mobility CFC Free PUF (Poly Urethane Foam) Insulation to ensure temperature stability and save power consumption Optional Accessories 7 Days circular temperature Recorder Multi Channel Data logging system Internal size (WxDxH) 58x46x46cms and outer size is 75x78x105 cms ( unit is also available with Data logger facility for maintaining the date in your PC on request)

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)