



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Portable Autoclave Single Drum Stainless steel

Item Code: S.T.R-06

The portable single drum Autoclave is highly efficient apparatus which is designed to meet the requirements of small dispensaries, nursing homes, laboratories, medical missions, clinics, ships, Operation theaters, Pharmaceuticals etc. It is very simple to operate and economical to use. The container is seamless made of thick stainless steel provided with safety valve, pressure gauge for the pressure inside the unit, steam releasing valve, electric heater running on 220V, 50Hz, the electric heater 3KW capacity, bakelite handle to carry the unit, unit based on heavy metal three leg stand, supplies with power cable and one spare gasket and one spare heating element.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com