



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Pulse Oximeter Table Model

Item Code:

We offer pulse oximeter of model, which is an extremely useful and comprehensive anesthetic monitoring device. It is widely used in family, hospital, oxygen bar, social medical organization, physical care in sports and others. Features physical and mechanical this is a portability, compact and desktop pulse oximeter measuring 22.0 cm h x 17.6 cm w x 7.5 cm d and weighing 720g without cable and accessories. Display the monitor screen is a lcd display that shows numeric patient information as well as alphanumeric status conditions. It can display the measurement results, instruction, troubleshooting messages, signal strength bar and time, its display parameters are %spo2 and pulse rate list of components 1 main unit, model 1 digital spo2 sensor for adult1 operator's manual 1 ac cable

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com