



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Revolving Stool

Item Code: HF-21

Revolving stool made of M.S frame duly chrome plated, the stool with five spokes top made of Good quality high density foam seat duly covered best quality of fabric/Leather type rexine stool has the provision of hydraulic height adjustment with the help of the lever which allow up and down of the seat position. units moves on the wheel.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)