



# HOSPITAL & MEDICAL EQUIPMENT



## Slit Lamp

Item Code:

Slit Lamp Max features are computer designed, multi-coated optics that provide unsurpassed viewing of ocular structures, the operator's field of view is maximized for precise visualization of ocular defects at the periphery.

### Specifications

#### Binocular Microscope

Eye Pieces : 10x, 15x  
Objective : 1x, 1.6x  
Total Magnification : 10x, 16x and 24x

#### Illumination Unit

Slit Image Rotation : 0 to 180°  
Tilting Illumination : Up to 20°  
Filter Disc : Cobalt Filter, Green Filter, Yellow Filter  
Natural Density and Open Aperture  
Slit Disphragm Disc : Wedge Shaped Diaphragm Of Infinitely  
Variable Slit Lengths

Halogen Lamp : 6 Volt 20 Watts

#### Standard Accessories

15x Eyepieces : Plastic Breath Shield  
Replacing Mirror : Hruby Lens  
Dust Cover and: Testing Rod

Fuse

Replacing Lamp one piece for illumination unit

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)