



Slit Lamp

Item Code:

Slit Lamp Max features are computer designed, multi-coated optics that provide unsurpassed viewing of ocular structures, the operator's field of view is maximized for precise visualization of ocular defects at the periphery.

Specifications

Binocular Microscope

Eye Pieces : 10x, 15x
Objective : 1x, 1.6x
Total Magnification : 10x, 16x and 24x

Illumination Unit

Slit Image Rotation : 0 to 180°
Tilting Illumination : Up to 20°
Filter Disc : Cobalt Filter, Green Filter, Yellow Filter
Natural Density and Open Aperture
Slit Disphragm Disc : Wedge Shaped Diaphragm Of Infinitely
Variable Slit Lengths

Halogen Lamp : 6 Volt 20 Watts

Standard Accessories

15x Eyepieces : Plastic Breath Shield
Replacing Mirror : Hruby Lens
Dust Cover and: Testing Rod

Fuse

Replacing Lamp one piece for illumination unit

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com