



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Stretcher Trolley Stainless steel with side rails

**Item Code:**

Stretcher trolley frame made of stainless steel round or square pipe head side adjustable with help of the lever provided on the head side, foot side fix provided with two handle to push the stretcher trolley top of the stretcher trolley foam padded duly covered with semi leather type material which is easy to clean and seepage proof, the unit moves on the wheels two are with break and two are without breaks. rubber shoe provided for stability on floor I.V. rod provision on the both side of the trolley. Side safety rails on both side

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)