



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Stretcher Trolley

Item Code:

Patient Stretcher Trolley with overall size 80"x22"x32" steel tubular frame work of trolley mounted on 15 cms dia swiveling castors wheels removable stretcher top made of S.S. sheet supported on Tubular for the oxygen gas cylinders and saline rod (Optional) finished with epoxy powder coated paint or chrome plating optional.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)