



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Surgical Cautery Machines (Digital)

Item Code: EM-01

surgical cautery of different model from 100 Watt to 400 watts, we offer best quality range of surgical cautery to our valuable clients. These are mainly used in surgeries of various kinds for the safety and security by using our machines. It is equipped with PEEC (plate electrode electronic control), which saves the patient from various dangers of burns, due to plate contact lack of age during surgery. These machines are very useful for the middle bipolar and monopolar electro-surgery. Input power required for operating it is 180-260 volts with maximum current of 3.5 amperes for 200 watt model. Both the bipolar and monopolar require high frequency leakage in 50ma rms and 100ma rms respectively. The power readouts work by actual power of +15% or 5 watts rated load. Indicators tones are adjustable to level of 45 db at 1 meter, whereas alarm tones are not adjustable and are set at 65 db minimum at 1 meter. The above specification is for 200 watt model. Detail specification of our other model is also on request. Model available up to 400 watt which also has feature of underwater cutting.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com