



An ISO 9001 & 14001 Co.

## HOSPITAL & MEDICAL EQUIPMENT



### Tonometer Schiotz

**Item Code:**

Rajas manufacture and supply first rate schiotz tonometer for the appropriate measurement of intraocular pressure. Our ophthalmic schiotz tonometer is highly apposite for complex and critical eye surgery. Our schiotz tonometer is available in superior quality brass parts and is also chrome plated. It has a red needle attached to it for the clear reading of the measured scale. We also make them available with conversion table in blue presentation case. So the schiotz tonometer for ophthalmologist is the most appropriate surgical equipment for the serious medical case related to the eyes. In addition to this, we are also one of the largest schiotz tonometer suppliers in India.

Tonometer according to Prof. Schiotz for measuring the Intraocular pressure. All main parts are made of High Quality Brass and are chrome plated. Red needle enables easy reading of the scale. The instrument comes complete with 3 weights (5.5 gram, 7.5gram, 10 gram) and conversion table in Presentation case.

## RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar  
Ambala Cantt- 133004 (Haryana) India  
E-mail : [rajasmedical@gmail.com](mailto:rajasmedical@gmail.com)  
[www.rajasindia.com](http://www.rajasindia.com)