



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Vaccine Refrigerator

Item Code: BB-16

Vaccine Refrigerator Capacity approx Vertical type with glass door. Temperature $^{\circ}\text{C} \pm 0,5\text{C}$. Accessories: Monitor for temperature with alarm visual and sound for high/low temperature, tropicalized. Auto defrosting. Adjustable shelves 2-3 in numbers, Lockable. 100% HCFC/CFC free. 240V/ 50 - 60Hz. 2-3. Inner made of Stainless steel 304 grade and outer made of M.S. sheet of heavy gauge duly powder coated paint finish, unit moves on the wheel to move around. Chart Recorder and Data Logger with RS 232 is Optional available on request. available in various capacity as per requirement.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com