



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Vertical Sterilizer

Item Code: S.T.R-03

we offer Vertical Sterilizer is an appliance to sterilize or clean the surfaces of the medical appliances from the bacteria, chemical deposits or dust particles. This appliance cleans the hospital and testing tools using high pressure steam. A sterilizer cleans testing jars, bottles, beakers, surgical instruments and hospital utensils from impurities and contaminations. We fabricate our entire gamut with quality approved materials.

Technical Specifications:

| | |
|----------------------|---|
| Construction | Sterilizer has inner and outer wall made of 304 grade stainless steel, middle jacket made of stainless steel or brass and gunmetal. |
| Technical Data | Steam Working Pressure 1.26 kgf/cm ² (2.2 kgf/cm ² in case High speed sterilizer) |
| Working Temperature | 121oC |
| Operating Voltage | 240 volts single phase, AC supply, 50 Hz. |
| Sterilization Period | 45 to 50 minutes |
| Components | Radiant locking system, cold water drainage system |
| Steam Exhaust | 5 to 7 provided with timer from 1 to 60 Min with double safety valve, steam releasing valve, pressure gauge for the pressure inside the chamber, and one pressure gauge to show the pressure of outer chamber both the chamber made of stainless steel lid of the unit made of heavy duty gun metal piece for extra safety, gasket, water level indicator to see the water inside, heavy duty power plug with socket provided with the unit |
| Accessories | one spare element, and two gasket , our Regular sizes are 900X500, and 1050X750. Other sizes also available on request. |

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com