



An ISO 9001 & 14001 Co.

HOSPITAL & MEDICAL EQUIPMENT



Wheel Chair Fixed

Item Code: HF-19

Over-all size Approx., 42"L x 28"W x 35"H, 18" Cushioned two Seat Folding & Back Fixed Detachable Armrest Front with 5" Solid Castors & Rear 26" Solid Wheels with Brakes. Arm rest cushion padded or PU made plastic handle to drive chair, foot rest made of plastic/ aluminum folding type, back seat has utility pocket, manual breaks provided. In the chair easy to fold and unfold the chair.

RAJAS ENTERPRISES (INDIA)

Nanhera Chowk, Nanhera, Post Office Kuldeep Nagar
Ambala Cantt- 133004 (Haryana) India
E-mail : rajasmedical@gmail.com
www.rajasindia.com